

Dépister le VIH en moins de 5 minutes maintenant  
c'est réalisable!

## RAPIDE

- Résultat en 60 secondes grâce à la technologie de l'ELISA rapide sur membrane

## PRATIQUE

- *Seulement 50 µl de sang total (capillaire ou Veineux EDTA), de sérum ou de plasma*

## SIMPLE

- Flacons unitaires pour chaque test
- Pas de minuteur
- Pas de compte-goutte
- Stockage à température ambiante
- Interprétation très simple (points)

## COMPLET

- Accessoires inclus (lancettes, pipettes, tampon alcool)

## FIABLE

- Excellentes sensibilité : 99,6 % et spécificité : 99,3 %
- Contrôle IgG au sein de chaque test



A photograph of two men embracing. The man on the left has dark, wavy hair and a beard, wearing a light-colored t-shirt. The man on the right has short dark hair and a beard, wearing a dark t-shirt. They are both smiling and looking at each other.

**INSTI**®

## HIV-1/HIV-2 Antibody Test

Introducing The World's Fastest  
Point-of-Care Test

- + Easy to use, early detection<sup>1</sup>
- + Clear and accurate results
- + Enhances patient flow
- + Expedites linkage to care



Is 60 seconds quick enough?

[biolytical.com](http://biolytical.com) • [insti.com](http://insti.com)

# INSTI<sup>®</sup> HIV-1/HIV-2 Antibody Test

Early IgM detection <sup>1</sup>


IgG Capture Procedural Control

No timers required


Lancet, pipette & alcohol swab included

“The ease of use of INSTI<sup>®</sup> and almost instant results make this test an accurate option for settings where recommended laboratory is not feasible...” - Adams, CDC Study <sup>1</sup>

## Sample & pour



## Read results immediately!



## Product Information

Information Type	Product Details
Method	Flow-through
Time to results	As little as 60 seconds
Storage conditions	35.6–86 °F (2–30 °C)
Shelf life	15 months
Sample type	Fingerstick, venous whole blood, plasma or serum

## Accuracy <sup>2</sup>

Sample Type	Sensitivity	Specificity
Venous Whole Blood	99.4%	99.7%
Plasma	99.5%	99.96%
Fingerstick Blood	99.6%	99.3%
Serum	99%	100%

## Ordering Information

### Product

INSTI<sup>®</sup> HIV-1/HIV-2 Antibody Test

INSTI<sup>®</sup> HIV-1 Positive/HIV-2 Positive HIV Negative Controls (3 x 1.0 mL vials)



### Trusted Technology That Works.

**Confidence.** The protein-A control feature serves as a built-in internal control and confirms correct sample addition, test procedure and proper test performance.

For ordering information, please contact:

1 604 204 6784  
info@biolytical.com

## References

<sup>1</sup> Adams S, et al. Performance evaluation of the point-of-care INSTI<sup>®</sup> HIV-1/2 Antibody Test in early and established HIV infections. *Journal of Clinical Virology* (2017) 91: 90–94.

<sup>2</sup> INSTI HIV-1/HIV-2 Antibody Test Package Insert (CAN), 51-1028 (J).